





CELLular
Dynamics
international




Humanizing the Tissue Chip: Use of Stem Cell Derived Tissues to Develop Biological Platforms

Kyle Kolaja, Ph.D., DABT, Fellow, ATS

July 2014



The Potential of iPS Cells: Genetic Diversity



Primary isolation of organs/tissues/cells

- >100 years since Harrison first cultured frog neurons
- >60 years since Gey first immortalized human cell line (HeLa)
- Immeasurable innovations, advancements, and knowledge

Yet, cell culture limitations haven't changed much and prevented the ultimate potential of replacing animal and human experiments

- Variability of isolation, timing, etc
- Degeneration of phenotype with time

IPS cell derived tissues have a number of advantages & improvements

- Footprint free method
- Human
- Gene editing/engineering
- Made from anyone



Primary Human Cells



Transformed Cell Lines



Primary cultured cells dedifferentiate and/or readily lose their phenotype in culture

- Primary hepatocytes
- Primary cardiomyocytes

Main driver for 3D culture models

- Ischemia-reperfusion stress induced during the isolation process
- Disruption of the tissue architecture and surviving in the new one

Maturation of stem cell derived tissues occurs in vitro


- Recapitulate embryonic development
- Incubate extended periods
- Electrical stimulation (cardiomyocytes)
- Small molecules
- Microphysiological systems

Make stem cell derived tissues mature faster!

- “fit for purpose”



- Largest Producer of hiPS and Derived Products
- Madison, WI, site Novato, CA (NASDAQ, ICEL)
- 150 employees (>600 FTE yrs of stem cell expertise)
- >800 Patents (Owned or Licensed)



● MesoDerm


- iCell Cardiomyocytes
- iCell Skeletal Myoblasts
- iCell Cardiac Progenitor
- iCell Hematopoietic Progenitor

● EctoDerm

- iCell Neurons
- iCell Dopa Neurons
- iCell Astrocytes

● EndoDerm

- iCell Hepatocytes
- iCell Endothelial Cells



● Custom iPS Cells (MyCell)



- **Robust manufacturing = enterprise wide quality management system**
 - Defined media and control of components
 - Substrate shift from feeder layers to recombinant proteins (e.g., laminin, vitronectin, etc)
 - Control of reagents from start to finish
 - Automation
- **Successful, broadly used items become commercialized → ACCESS**
 - Media and substrates above
 - Micro-arrays great example
 - Academia – Govt – Industry
 - Homemade to QC product

Journal of Biomedical Discovery and Collaboration



Case Study

The emergence and diffusion of DNA microarray technology

Tim Lenoir**† and Eric Giannella†

Open Access

Address: Jackson Laboratory for New Technologies in Society, Duke University, John Hope Franklin Center, 2204 Erwin Road, Durham, North Carolina 27708-0402, USA


Email: Tim Lenoir* - lenoir@duke.edu; Eric Giannella - ericgiannella@duke.edu

* Corresponding author † Equal contributors



Scale-Up Manufacturing

- Quality
- Quantity
- Purity



CDI Manufacturing Benchmarks

(cells per day, >95% purity)

- 2 billion iPS cells
- 1 billion iCell Cardiomyocytes
- 1 billion iCell Neurons
- 0.5 billion iCell Endothelial cells
- 0.4 billion iCell Hepatocytes
- ...

Scale-Out Manufacturing

- 1000's of individuals
- Billions of cells

NHLBI Next Generation Genetic Association Studies (RFA-HL-11-066)

- 250 patient samples - HyperGEN cohort
- GWAS – Left Ventricular Hypertrophy (LVH)
- Derive iPS cells and cardiomyocytes from all 250 individuals
- Induce hypertrophy phenotype, perform molecular analyses
- Correlate GWAS findings with in vitro phenotype



- **California Institute for Regenerative Medicine (CIRM)**
- **Human iPS Cell Initiative – 3 Awards (Total \$32M)**
 - Sample Collection (7 awardees - \$0.5M - \$1.5M)
 - iPS Cell Derivation (CDI - \$16M)
 - iPS Cell Banking (Coriell - \$10M; CDI primary subcontractor)
- **iPS Cell Derivation**
 - 3000 donors (healthy & disease phenotypes)
 - 3 iPS cell clones per donor
 - Disease categories: epilepsy, autism, cerebral palsy, cardiomyopathy, Alzheimer's disease, eye diseases, hepatitis (HCV), non-alcoholic steatohepatitis (NASH), pulmonary fibrosis
 - Derived from peripheral blood (preferred) or skin fibroblasts
 - Episomal "footprint-free" method
- **CDI – Coriell Partnership**
 - Brings together expertise in electronic record-keeping, sample tracking, iPS cell derivation & characterization, cell banking & distribution
 - Joint facility located within the Buck Institute of Aging, Novato, CA


CALIFORNIA INSTITUTE FOR
REGENERATIVE MEDICINE




CELLular
Dynamics
international




CORIELL INSTITUTE
FOR MEDICAL RESEARCH







- **First stem cell derived cell type**
 - ~4-5 years of publications now
 - Used in regulatory filings to support claims
- **Improvements in purity and quantity sped the growth and adoption**
- **Proof of comparability (+) established**
 - Gene expression
 - Morphology
 - Electrophysiology and contractility
 - Biochemical properties
 - Functional (pharm and tox)
- **Rapidly emerging opportunity in arrhythmia detection, but ample applications in pharmacology, toxicology, and disease biology research**






Gene Expression



Electrophysiology





IB20= lowest tested concentration resulting in 20% irregular beats



Human Cardiomyocyte Arrhythmia Assay Validation: 23 compounds with known *in vivo* effect

- 12 Pro-arrhythmic
- 11 Non-arrhythmic
- **IB20 30 μ M**
 - One False Positive
 - No False Negatives

Drug	IB20 [®] (μ M)	hERG	QT	Clinical [®] arrhythmia
Dofetilide	0.003	(+)	(+)	(+)
Ouabain	0.03	(-)	(-)	(+)
Aconitine	0.03	(-)	(-)	(+)
Cisapride	0.03	(+)	(+)	(+)
E-4031	0.03	(+)	(+)	(+)
Astemizole	0.03	(+)	(+)	(+)
Terfenadine	0.3	(+)	(+)	(+)
Flecainide	1	(+)	(+)	(+)
<i>Alfuzosin</i>	1	(-)	(+)	(-)
Thioridazine	3	(+)	(+)	(+)
Quinidine	10	(+)	(+)	(+)
Erythromycin	30	(+)	(+)	(+)
Sotalol	30	(+)	(+)	(+)
Fluoxetine	>30	(+)	(+)	(-)
Verapamil	>30	(+)	(\pm)	(-)
Moxifloxacin	>100	(+)	(+)	(+)
Amiodarone	>100	(+)	(+)	(+)
Ranolazine	>100	(+)	(+)	(-)
Captopril	>100	(-)	(-)	(-)
Rofecoxib	>100	(-)	(-)	(-)
Amoxicillin	>1000	(-)	(-)	(-)
Aspirin	>1000	(-)	(-)	(-)
Nifedipine	>3	(-)	(-)	(-)

- **IB20= lowest tested concentration resulting in 20% irregular beats**

Guo et al. Toxicol Sci. 2011 Sep;123(1):281-9.



Human Cardiomyocyte Arrhythmia (hCAR) Assay

2nd Set of Validation and Model Refinement

Drug	C ₅₀ (nM)	IB ₅₀ (µM)	hERG	QT	QTp	in vivo ECG			hERG			
Digitoxin	33	0.003	(-)	(-)	(+)	32	Cyclosporin A	1,458	100	(-)	(-)	(-)
Dofetilide	6	0.003	(+)	(+)	(+)	14	Lidocaine (i.v.)	36,000	100	(+)	(-)	(-)
Digoxin	3	0.01	(-)	(-)	(+)	3	Pimobendan	164	100	(-)	(-)	(-)
Ouabain	170	0.01	(-)	(-)	(+)	3	Ranolazine	6,005	100	(+)	(+)	(-)
Aconitine	77	0.03	(-)	(-)	(+)	3	Nifedipine	154	>30	(-)	(-)	(-)
Artemizole	8	0.03	(+)	(+)	(+)	3	Amiflopridine	793	>30	(-)	(+)	(-)
E-4031	15	0.03	(+)	(+)	(+)	3	Epridil	3,258	>30	(+)	(+)	(+)
Pimozide	217	0.1	(+)	(+)	(+)	3	Ceftriaxone	300	>30	(-)	(-)	(-)
Serinidole	318	0.1	(+)	(+)	(+)	0	Cibenzoline	2,168	>30	(+)	(+)	(-)
Cisapride	125	0.3	(+)	(+)	(+)	0	Calcitripridine	453	>30	(-)	(-)	(-)
Geldanamycin	16,300	0.3	(-)	(-)	(+)	0	Diltiazem	552	>30	(-)	(-)	(-)
Ixabucin	123	0.3	(-)	(+)	(+)	0	Difenhydramine	157	>30	(-)	(+)	(-)
Terfenadine	300	0.3	(+)	(+)	(+)	0	Fluorouracil	4,613	>30	(-)	(-)	(+)
Alfuzosin	56	1	(-)	(-)	(-)	0	Fluoxetine	425	>30	(+)	(+)	(-)
Dobutamine	3,319	1	(+)	(-)	(-)	0	Imipramine	1,070	>30	(+)	(+)	(-)
Doxorubicin	15,344	1	(-)	(+)	(+)	0	Ketoconazole	17,628	>30	(+)	(+)	(-)
Flecainide	1,531	1	(+)	(+)	(+)	10	Loratadine	23	>30	(-)	(-)	(-)
Pentamidine	2,181	1	(-)	(+)	(+)	10	Nitrendipine	150	>30	(-)	(-)	(-)
Tacrine	100	1	(-)	(-)	(-)	10	Olanzapine	74	>30	(-)	(-)	(-)
Amphotericin B	33,318	3	(-)	(+)	(+)	10	Rosiglitazone	1,673	>30	(-)	(-)	(-)
Arsenic Trisulfide	12,132	3	(-)	(+)	(+)	10	Troglitazone	6,387	>30	(-)	(-)	(-)
Clozapine	1	3	(-)	(+)	(-)	10	Verapamil	315	>30	(+)	(-)	(-)
Mitomycin	3,311	3	(-)	(+)	(+)	10	Acepromidophenol	130,000	>100	(-)	(-)	(-)
Prentalamine	70	3	(-)	(+)	(+)	10	Aldidem	284	>100	(-)	(-)	(-)
Sunitinib	253	3	(+)	(+)	(+)	10	Amlodipine	3,374	>100	(+)	(+)	(-)
Thioridazine	1,731	3	(+)	(+)	(+)	31	Atenolol	1,284	>100	(-)	(-)	(-)
Zimelidine	328	3	(+)	(-)	(-)	10	Captopril	2,466	>100	(-)	(-)	(-)
Ajmaline (i.v.)	105	10	(+)	(+)	(+)	10	Colchicine	16	>100	(-)	(-)	(-)
Chlorpromazine	2,630	10	(+)	(+)	(+)	10	Cyclophosphamide	153,200	>100	(-)	(-)	(-)
Carthromycin	6,029	10	(+)	(+)	(+)	10	Desazotane	136,052	>100	(-)	(-)	(-)
Carboline	7,500	10	(-)	(-)	(-)	10	Levomepromidol	136	>100	(-)	(-)	(-)
Desipramine	601	10	(+)	(+)	(+)	10	Mechlorethamine	?	>100	(-)	(-)	(-)
Eprubicin	16,036	10	(-)	(+)	(+)	10	Moxifloxacin	10,276	>100	(+)	(+)	(-)
Nefazodone	4,358	10	(+)	(+)	(-)	10	Nimetazide	15,000	>100	(-)	(-)	(-)
Phentolamine	100	10	(-)	(-)	(-)	10	Pemoline	6,525	>100	(-)	(-)	(-)
Quinidine	21,573	10	(+)	(+)	(+)	01	Rofecoxib	1,021	>100	(-)	(-)	(-)
Erythromycin (i.v.)	34,064	30	(+)	(+)	(+)	10	Tolcapone	21,555	>100	(-)	(-)	(-)
Fluvoxamine	1,257	30	(+)	(-)	(-)	10	Zalcitabine	115	>100	(-)	(-)	(-)
Imatinib	3,541	30	(-)	(+)	(-)	10	Amoxicillin	17,036	>1000	(-)	(-)	(-)
Mefenazine	11,161	30	(+)	(-)	(-)	30	Aspirin	10,000	>1000	(-)	(-)	(-)
Propafenone	4,327	30	(+)	(+)	(+)	00	-	-	-	4.9	-	-
Propranolol (i.v.)	183	30	(-)	(+)	(-)	00	-	-	-	30	-	-
Sotalol	14,733	30	(+)	(+)	(+)	00	-	-	-	4.0	-	-
			Z/	>100	>100	-	-	-	-	4.3	-	-
			28	>100	>100	-	-	-	-	30	-	-
			29	>100	>100	-	-	-	-	>300	-	-
			30	>100	>100	-	-	-	-		-	-

83 Compounds

~82% -- arrhy. prediction
>90% -- QT prediction


30 Internal Compounds

80% -- arrhy. prediction
95% -- QT prediction


Refining the Human iPSC-Cardiomyocyte Arrhythmic Risk Assessment Model

Liang Guo,¹ Luke Coyle, Rory M.C. Abrams,² Raymond Kempczak,² Eric T. Chiao,⁴ and Kyle L. Kolajczak²


- **Product launch → regulatory evaluation in 3 years**
- **iPS cell-derived cardiomyocytes are being evaluated for use in arrhythmia assessment & as a replacement for thorough QT studies**




Nature Reviews Drug Discovery (Aug, Sept 2013)




BoNT Receptor Protein Expression



BoNT Receptor Gene Expression



BoNT Receptor Cleavage



- iCell Neurons express the receptors and enzymatic targets necessary for BoNT cell entry and catalytic activity
- iCell Neurons reproducibly show equivalent or greater sensitivity to BoNT activity vs. rat spinal cord cells


- Assess the potency of botulinum neurotoxin (BoNT) better than rat spinal cord neurons or mouse LD50.
- a consortium of BoNT manufacturers is in the process of validating the **use of Stem Cell derived Neurons** to replace the current industry “gold” standard, a high-cost and labor-intensive in vivo bioassay.

Novel Application of Human Neurons Derived from Induced Pluripotent Stem Cells for Highly Sensitive Botulinum Neurotoxin Detection

Regina C. M. Whitmarsh,* Monica J. Strathman,† Lucas G. Chase,* Casey Stankewicz,† William H. Tepp,*
Eric A. Johnson,* and Sabine Pellett^{†,1}

*Department of Bacteriology, University of Wisconsin, Madison, Madison, Wisconsin 53706 and †Cell Biology Group, Cellular Dynamics International, Inc., Madison, Wisconsin 53711





iPSC technology can be used to model human **Innate**, **Induced** and **Infectious** diseases that cannot be interrogated using conventional cell lines, primary cells or animal models

Adapted from Grskovic, et al. (2011)


Reviews & summaries of disease-specific iPSCs created:

- Grskovic, et al. (2011) *Nature Reviews Drug Discovery*
- Rajamohan, et al. (2012) *Bioessays*
- Trounson, et al. (2012) *Current Opinion Genetics & Development*




Cardiac Hypertrophy

Mouse Phenotype (increased cell size)




Normal Disease



Endothelin Induced iCell CM





Control ET-1 (10 nM)




Flow Cytometry



qRT-PCR







Robust and Reproducible in-vitro assays (BNP-based readouts)




Steatohepatitis

Amiodarone Induced iCell HC




iCell Hepatocytes HCV Infection (Clinical Genotypes)

Luc Expressing HCV
pseudoparticle
(HCVpp) uptake




HCVpp encoding Firefly luciferase


Inhibition of HCVpp Uptake
by anti-CD81 Ab



**iCell Hepatocytes are Susceptible to Multiple
HCV Genotypes**




HCVcc - Cell Culture Passaged
Virus (Genotype 1a/2a)




Patient Serum HCV
(Genotype 1a)

iCell Neurons Physiologic VZV Latent Infection



VZV infection did not produce a cytopathic effect



Viral Particles and Capsids in iCell Neurons


Yu, et al. (2013) J Neurovirology



iCell Neurons: Novel Alzheimer's Disease Biology & Use in Phenotypic Screens

Novel Alzheimer's Disease Biology

Tau Secretion



◀ **Constitutive extracellular secretion of tau protein by an unconventional mechanism. Secretion was influence by both time and temperature.**




Percentage of tau released from cells.

	Tau in medium (pg/ml)	Intracellular tau (ng/ml)	% Release
HEK 293 Trex tau wt	134 ± 19	102 ± 18	0.1
HEK 293 Trex tau ΔK280	261 ± 36	267 ± 86	0.1
iCell® Neurons	59 ± 6	20 ± 2	0.3


Chai X, et al. (2012) Neurobiol Dis.

β-Amyloid Processing

CDI (2013)





▶ **Endogenous levels of Aβ40 was modulated upon treatment with a secretase inhibitor.**





Relevant Neuronal Model for Alzheimer's Disease




▼ **iCell Neurons showed sensitivity to Aβ1-42.**



◀ ◀ **Neuroprotection from induced Aβ1-42 toxicity was demonstrated by small-molecule inhibition and shRNA knockdown of CDK2.**




Scatter Plot



▶ **HT phenotypic screen identified 19 hits, including a CDK2 inhibitor confirming reliability and sensitivity of the platform.**

Xu X, et al. (2013) Stem Cell Res



Neuronal Diseases

Amyotrophic lateral sclerosis
Spinal muscular atrophy
Olivopontocerebellar atrophy
Parkinson's disease
Huntington's disease
Down's syndrome
Fragile X syndrome
Friedrichs Ataxia
Familial dysautonomia
Rett's syndrome
Mucopolysaccharidosis type IIIB
Schizophrenia
X-linked adrenoleukodystrophy
childhood cerebral ALD
Adrenomyeloneuropathy
Autism spectrum disorders
Angelman syndrome
Pradder-Willi

Skin


Recessive dystrophic
epidermolysisbullosa

Eye

Retinitis pigmentosa
Age-related cataract
Gyrate atrophy

Multi-organ

Down syndrome - Trisomy 21
Shwachman-Bodian-Diamond syndrome
Dyskeratosiscongenita



Current status of drug screening and disease modelling in human pluripotent stem cells

*Divya Rajamohan, Elena Matsa, Spandan Kalra, James Crutchley, Asha Patel, Vinoj George and Chris Denning**

Bioessays 35: 281-298, © 2012 WILEY Periodicals, Inc.

Muscle

Duchene Muscular Dystroph
Becker muscular dystrophy
Hutchinson-Gilford progeria syndrome

Metabolic

Gaucher disease type III
Lesch-Nyhan syndrome
Juvenile Diabetes
Type 2 diabetes
Familial hypercholesterolemia
Alpha1-antitrypsin deficiency
Glycogen storage disease type 1a

Immune

Adenosine deaminase deficiencyassociated
severe combined
immunodeficiency (ADA-SCID)
Multiple Sclerosis

Cardiovascular Diseases

Flavors of long QT syndrome
CPTV
LEOPARD syndrome
Timothy Syndrome

Haematological

Sickle cell anaemia b-Globin alleles
Fanconi anaemia
Acquired myeloproliferativedisordes
b-Thalassaemia major (Cooley's
anaemia)



Historical Phenotypic Screens

- MMOA knowledge not required
- Target ID difficult
- Lack of relevant human biology


1980-2000 Target-based Screens

- Driven by genomics & informatics
- Limited new validated targets under study
- Lack of translatability in vivo
- **↑ R&D investments; ↓ first-in-class drugs**

2000-Present Phenotypic Screens

- Rapid discovery of disease modifying molecules
- Improved Target ID (genetic and chemical proteomics)
- **Improved relevance of human biology**

Swinney & Anthony (2011) Nat Rev Drug Discovery



ANALYSIS

How were new medicines discovered?

David C. Swinney^{1*} and Jason Anthony^{2*}

Abstract | Preclinical strategies that are used to identify potential drug candidates include target-based screening, phenotypic screening, modification of natural substances and biologic-based approaches. To investigate whether some strategies have been more successful than others in the discovery of new drugs, we analysed the discovery strategies and the molecular biology that were used to identify potential drug candidates in 2008. Out of the 100 most successful strategies, 17 were biologic-based approaches, and 83 were target-based approaches. The results also show that 17 of these strategies were the major focus for first-in-class drug discovery, and 17 of these strategies contribute to the research and development of new drugs.

SciBX
Science-Business eXchange


Phenotypic screening, take two
By Joanne Kotz, Senior Editor

Kotz (2012) SciBx – April 12, 2012



iPSC technology can be used to model diseases with known MMOA, as well as in phenotypic-based screens for complex diseases with unknown genetic mechanisms





Regenerative Medicine and the Future of iPSC: Cell, Tissue and Organ Creation



Trachea Brain



Patent pending, Hideki Taniguchi, Takanori Takebe: PCT/JP2012/074840



IPS-derived cells improvement over primary culture

- Amenable to genetic engineering
- Maturing phenotype
- Relevant disease models can be induced or derived

Improved functionality → Ask better questions

Robust manufacturing a necessity

iPS cells allow direct control over genetic diversity

- Patient disease phenotype recapitulation in vitro
 - Retrospective clinical trials
 - Prospective clinical trials ?
- **Clinical applications have potential to completely change medicine**

