

---

**Lentiviral vectors and protocols for creation of stable hESC lines for fluorescent tracking and drug resistance selection of cardiomyocytes.**

**Journal:** PLoS One

**Publication Year:** 2009

**Authors:** Hiroko Kita-Matsuo, Maria Barcova, Natalie Prigozhina, Nathan Salomonis, Karen Wei, Jeffrey G Jacot, Brandon Nelson, Sean Spiering, Rene Haverslag, Changsung Kim, Maria Talantova, Ruchi Bajpai, Diego Calzolari, Alexey Terskikh, Andrew D McCulloch, Jeffrey H Price, Bruce R Conklin, H S Vincent Chen, Mark Mercola

**PubMed link:** 19352491

**Funding Grants:** Chemical Genetic Approach to Production of hESC-derived Cardiomyocytes, Development of Neuro-Coupled Human Embryonic Stem Cell-Derived Cardiac Pacemaker Cells., Burnham Institute CIRM Stem Cell Training Grant (Type II)

**Public Summary:**

**Scientific Abstract:**

**BACKGROUND:** Developmental, physiological and tissue engineering studies critical to the development of successful myocardial regeneration therapies require new ways to effectively visualize and isolate large numbers of fluorescently labeled, functional cardiomyocytes. **METHODOLOGY/PRINCIPAL FINDINGS:** Here we describe methods for the clonal expansion of engineered hESCs and make available a suite of lentiviral vectors for that combine Blasticidin, Neomycin and Puromycin resistance based drug selection of pure populations of stem cells and cardiomyocytes with ubiquitous or lineage-specific promoters that direct expression of fluorescent proteins to visualize and track cardiomyocytes and their progenitors. The phospho-glycerate kinase (PGK) promoter was used to ubiquitously direct expression of histone-2B fused eGFP and mCherry proteins to the nucleus to monitor DNA content and enable tracking of cell migration and lineage. Vectors with T/Brachyury and alpha-myosin heavy chain (alphaMHC) promoters targeted fluorescent or drug-resistance proteins to early mesoderm and cardiomyocytes. The drug selection protocol yielded 96% pure cardiomyocytes that could be cultured for over 4 months. Puromycin-selected cardiomyocytes exhibited a gene expression profile similar to that of adult human cardiomyocytes and generated force and action potentials consistent with normal fetal cardiomyocytes, documenting these parameters in hESC-derived cardiomyocytes and validating that the selected cells retained normal differentiation and function. **CONCLUSION/SIGNIFICANCE:** The protocols, vectors and gene expression data comprise tools to enhance cardiomyocyte production for large-scale applications.

---

**Source URL:** <http://www.cirm.ca.gov/about-cirm/publications/lentiviral-vectors-and-protocols-creation-stable-hesc-lines-fluorescent>