
Human iPSC derived progenitors bioengineered into liver organoids using an inverted colloidal crystal poly (ethylene glycol) scaffold.

Journal: Biomaterials

Publication Year: 2018

Authors: Soon Seng Ng, Kouros Saeb-Parsy, Samuel J I Blackford, Joe M Segal, Maria Paola Serra, Marta Horcas-Lopez, Da Yoon No, Sotiris Mastoridis, Wayel Jassem, Curtis W Frank, Nam Joon Cho, Hiromitsu Nakauchi, Jeffrey S Glenn, S Tamir Rashid

PubMed link: 30149262

Funding Grants: Generation of functional cells and organs from iPSCs

Public Summary:

Generation of "mini" organs-in-a-dish (organoids) from human pluripotent stem cells represents a powerful new approach to study organ development and model disease. Here, we have developed a full-defined protocol to generate liver organoids, and use this system to identify signaling pathways involve in liver development.

Scientific Abstract:

Generation of human organoids from induced pluripotent stem cells (iPSCs) offers exciting possibilities for developmental biology, disease modelling and cell therapy. Significant advances towards those goals have been hampered by dependence on animal derived matrices (e.g. Matrigel), immortalized cell lines and resultant structures that are difficult to control or scale. To address these challenges, we aimed to develop a fully defined liver organoid platform using inverted colloid crystal (ICC) whose 3-dimensional mechanical properties could be engineered to recapitulate the extracellular niche sensed by hepatic progenitors during human development. iPSC derived hepatic progenitors (IH) formed organoids most optimally in ICC scaffolds constructed with 140µm diameter pores coated with type I collagen in a two-step process mimicking liver bud formation. The resultant organoids were closer to adult tissue, compared to 2D and 3D controls, with respect to morphology, gene expression, protein secretion, drug metabolism and viral infection and could integrate, vascularise and function following implantation into livers of immune-deficient mice. Preliminary interrogation of the underpinning mechanisms highlighted the importance of TGFbeta and hedgehog signalling pathways. The combination of functional relevance with tuneable mechanical properties leads us to propose this bioengineered platform to be ideally suited for a range of future mechanistic and clinical organoid related applications.

Source URL: <https://www.cirm.ca.gov/about-cirm/publications/human-ips-derived-progenitors-bioengineered-liver-organoids-using-inverted>